



Examination lamps

Application

- * ENT,gynecology and obstetrics department , outpatient department examination;
- *Used as assisted illumination in operating room.



Technical parameters

Model	YJ-JC01
Bulb	Halogen
Supply Voltage	220V
Power	35W
Color temperature	4000K
Illumination(50cm)	30000Lux
Spot size adjust	×
Brightness adjust	×
Height adjust	√
Lifetime	4000h
Packing size	93*48*10cm
N.W/G.W	4KG/5KG



Examination lamp

Product Features

- *Mini lamp holder with diameter of 75mm;
- *Illumination using Osram imported from Germany, with stable light and long working time;
- *Flexible arm can freely position,easy to move;
- *Even bright white light.



Technical parameters

Model	YJ-JC02
Bulb	LED
Supply Voltage	220V
Power	3W
Color temperature	4500K
Illumination(50cm)	16000Lux
Spot size adjust	×
Brightness adjust	×
Height adjust	√
Lifetime	50000h
Packing size	93*48*10cm
N.W/G.W	4KG/5KG



Examination lamp

Product Features

- *Mini lamp holder with diameter of 75mm;
- *Illumination using Osram imported from Germany, with stable light and long working time;
- *Flexible arm can freely position,easy to move;
- *Even bright white light.



Technical parameters

Model	YJ-JC02 (With battery)
Bulb	LED
Supply Voltage	220V
Power	3W
Color temperature	4500K
Illumination(50cm)	16000Lux
Spot size adjust	×
Brightness adjust	×
Height adjust	√
Lifetime	50000h
Packing size	93*48*10cm
N.W/G.W	4KG/5KG



Examination lamp

Product Features

- * Single-hole cold light, vertical and moveable.
- * Adjustable height between 1000- 1600mm.
- * With the spherical joint, the lamp holder is upright and side-ling adjustable.
- * Pedal switch is installed in the base, which makes the lamp has more convenient use.

Application

- * ENT, Urology, Obstetrics and Gynecology;
- * Used as assisted light in the operating theatre.

Technical parameters

Model	YJ-JC05
Bulb	Halogen
Supply Voltage	220V
Power	25W
Color temperature	3000 ± 500K
Illumination(50cm)	≥ 12000Lux
Spot size adjust	×
Brightness adjust	×
Height adjust	√
Lifetime	4000h
Packing size	115*28*24cm
N.W/G.W	10KG/11KG



Examination lamp

Product Features

- * It adopts imported LED light source, which has high brightness, long life and better color rendering.
- * Adjustable height between 1000- 1600mm.
- * The lamp head can be adjusted 90° front and rear with a ball joint.
- * Brightness adjustment is available as needed.
- * The foot switch is mounted on the base and is easy to use

Application

- * ENT, Urology, Obstetrics and Gynecology;
- * Used as assisted light in the operating theatre.

Technical parameters

Model	YJ-JC06
Bulb	LED
Supply Voltage	220V
Power	8*3W
Color temperature	5000 ± 200K
Illumination(50cm)	≥80000Lux
Spot size adjust	×
Brightness adjust	√
Height adjust	√
Lifetime	50000h
Packing size	116*28.5*24cm
N.W/G.W	10KG/11KG



Examination lamp

JSF-JC series examination lamps are designed for the use in outpatient examination and small surgery.

Product Features

- *Mini lamp holder with diameter of 75mm;
- *Illumination using Osram imported from Germany, with stable light and long working time;
- *Flexible arm can freely position, easy to move;
- *Even bright white light.

Application

- * ENT, gynecology and obstetrics department, outpatient department examination;
- * Used as assisted illumination in operating room.

Technical parameters

Model	YJ-JCA02
Bulb	LED
Supply Voltage	AC110-240V
Power	9W
Color temperature	6500K
Illumination(50cm)	40000Lux
Spot size adjust	×
Brightness adjust	√
Height adjust	√
Lifetime	50000h
Packing size	93*48*10cm
N.W/G.W	4.5KG/5.5KG

Examination lamp



JSF-JC series examination lamps are designed for the use in outpatient examination and small surgery.

Product Features

- *Mini lamp holder with diameter of 75mm;
- *Illumination using Osram imported from Germany, with stable light and long working time;
- *Flexible arm can freely position, easy to move;
- *Even bright white light.

Application

- * ENT, gynecology and obstetrics department, outpatient department examination;
- * Used as assisted illumination in operating room.

Technical parameters

Model	YJ-JCA02
Bulb	LED
Supply Voltage	AC110-240V
Power	9W
Color temperature	6500K
Illumination(50cm)	40000Lux
Spot size adjust	×
Brightness adjust	√
Height adjust	√
Lifetime	50000h
Packing size	41.5*31.5*12.5cm
N.W/G.W	3KG/4KG



Examination lamp

Three high quality LED bulbs
Illumination 40000Lux



Technical parameters

Model	YJ-JCU02
Bulb	LED
Supply Voltage	AC110-240V
Power	9W
Color temperature	6500K
Illumination(50cm)	40000Lux
Spot size adjust	×
Brightness adjust	×
Height adjust	√
Lifetime	50000h
Packing size	93*48*10cm
N.W/G.W	4KG/5KG



Examination lamp

Product Features

- * High quality LED modules;
- * Spot size and brightness can be adjusted based on different needs;
- * The high-quality metal hose can be bent freely, making it easier and more flexible to use.

Application

- * Outpatient examinations of ENT, obstetrics and gynecology, dentistry, etc;
- * Local lighting for minor operations;
- * Auxiliary lighting in operating room.

Technical parameters

Model	YJ-JCA03
Bulb	LED
Supply Voltage	AC110-240V
Power	12W
Color temperature	5000K
Illumination(50cm)	9000Lux
Spot size adjust	√
Brightness adjust	√
Height adjust	×
Lifetime	50000h
Packing size	57*47*18cm
N.W/G.W	4.5KG/5.5KG



Examination lamp

Product Features

- * High quality LED modules;
- * Spot size can be adjusted based on different needs;
- * The high-quality metal hose can be bent freely, making it easier and more flexible to use.
- * Foot switch, easy to use.



Foot Switch


Technical parameters

Model	YJ-JCU03
Bulb	LED
Supply Voltage	AC110-240V
Power	12W
Color temperature	5000K
Illumination(50cm)	9000Lux
Spot size adjust	√
Brightness adjust	×
Height adjust	×
Lifetime	50000h
Packing size	57*47*18cm
N.W/G.W	4KG/5KG



Examination lamp

This type of examination lamp mainly Applies to gynecology and obstetrics Examination surgery, using LED lighting, the spot t boundary is very clear and the light is even, it uses the guide rail to clip on obstetric table, hose can be bent easily, has very convenient use.



Technical parameters

Model	YJ-JCA04
Bulb	LED
Supply Voltage	AC110-240V
Power	5W
Color temperature	5000K
Illumination(50cm)	10000Lux
Spot size adjust	×
Brightness adjust	×
Height adjust	×
Lifetime	20000h
Packing size	57.5*37*9.5cm
N.W/G.W	3KG/4KG