



Henan Forever Medical CO.,LTD www.yongjiumed.com

YJ-ECG12 Digital Electrocardiograph Series

Quality assurance/safe and effective/customer service

YJ-ECG12 series

digital electrocardiograph



YJ-ECG12 digital electrocardiograph

Features

10.1" color LCD screen, color text menu, man-machine interface, which is convenient and intuitive 12 leads synchronous acquisition, 12 ECG

Waveforms display in same screen and synchronous printing

Touch screen operation, which makes the operation more quick and simple, time and labor saving

Manually and automatically print the standard 12 leads, 3 leads plus rhythm leads, 6 leads, 6 leads plus rhythm leads and 60s arrhythmia analysis. Large volume ECG data storage, and capacity expansion of SD storage card supported

Electrode loss alarm and indication of electrode loss parts

[Email:forermed@yjcompany.com](mailto:forermed@yjcompany.com)

[Tel:+86-37161315851](tel:+86-37161315851)



The V1, V3 and V5 chest lead examination mode is easy to operate, which is applicable for the large scale physical examination.

It takes the unique high precision digital filter to eliminate the base line drift and other interferences, and the ECG waveform doesn't distort, to improve the base line drift capacity greatly, and facilitate the judgment of the waveform.

It can store 250 ECG data files

Data storage (optional)



Large volume ECG data storage, capacity expansion supported

Interface



Color text menu, man-machine interface, which is convenient and intuitive

Fittings



Lead wire



clamp



suction bulb

Touch screer



Quick and simple to operate, time and labor saving

The printed picture is only for reference for there may be slight color deviation between it and the product. The color of the physical products shall prevail.