

YJ-V12

Veterinary multi-parameter monitor



Specially designed for pets hospital to measure accurately, simple and friendly operating display interface. Meet most of the needs of pets hospital operating room

Features

- 1. 12.1 inch high brightness color LCD/touch screen;
- 2. Various operation modes of buttons, knobs and touch screens;
- 3. It can be used in cats, dogs, pigs, horses and other small, medium and large animals.
- 4. Support multi-channel 8 waveform with screen display;
- 5. Voice, text and visual alarm function, level 3 alarm
- 6. Multiple interface display: standard interface, trend coexistence interface, oxyCRG interface, large font interface;
- 7. Supports multiple languages: Chinese, English, Spanish, polish, Turkish, French, Portuguese, Russian;
- 8. Built-in high capacity battery, long battery life, easy to transfer animals;
- 9. Optional print output function, support alarm trigger print;
- 10. Anti - defibrillation, anti - frequency electric knife interference. 11. Support for wired or wireless central monitoring systems;
- 12. Optional extension function module (dual channel IBP, ETCO2), plug and play, software automatic recognition, interface dynamic adjustment, The module can be hot-plug;
- 13. Animal special accessories for different types of animals;



Attachment

- 12.1 inch high brightness color LCD/touch screen;
- Built-in high capacity battery
- Standard configuration:
SPO2, ECG, TEMP, PR/HR, NIBP, RESP
- Options:
ETCO2, IBP, Nellcor SPO2, Suntech NIBP, recorder, trolley
- Six functions



ECG lead wire



Blood oxygen probe



Blood pressure cuff



Blood oxygen extension cord



Blood pressure extension cord



Ground wire

Specially designed for pets hospital o

- Reliable and comprehensive multi-parameter monitor
- Blood pressure extension cord
- Ground wire
- Probe, Including SPO2, ECG, TEMP, PR/HR, NIBP, RESP
- Exclusive patented BMV blood pressure measurement technology.
- Anti - defibrillation, anti - frequency electric knife interference.
- Support for wired or wireless central monitoring systems;
- System software supports network upgrade



standard interfac



oxyCRG interface



large font interface



trend coexistence interface

